

Abstract

Objectives. Sexually transmitted infections (STI) are widely abundant all over the world, presenting constantly a yearly growing trend and consequently increased annual health costs. STIs may cause acute symptoms, chronic infections, infertility, ectopic pregnancy, blindness, cervical cancer, susceptibility to HIV infection and death. Most cases are treatable, however no or late treatment may lead to severe health complications. This is further emphasized in view of the common carriage of STIs in asymptomatic state. Therefore, there is a need for detection of early stages of STIs in order to control and prevent these diseases. Savyon Diagnostics has recently been engaged with development of a novel molecular-based screening test for simultaneous detection of the most abundant STIs utilizing its proprietary NanoCHIP®XL molecular electronic microarray system. The panel of pathogens is composed of *Chlamydia trachomatis*, *Neisseria gonorrhoeae*, *Trichomonas vaginalis*, *Ureaplasma urealyticum/parvum*, *Mycoplasma genitalium* and *Mycoplasma hominis*. The *Ureaplasma urealyticum/parvum* results are reported in a semi-quantitative manner, in order to better differentiate between colonization and acute states, and by that avoiding unnecessary antibiotic treatment. The aim of this work is to demonstrate the utility of the newly developed test for screening purposes of STIs especially in high-risk populations, demonstrating a highly performing, time-saving and cost-effective test.

Methods. DNA was extracted from characterized urethral samples or urine swab specimens using a variety of readily available manual and automatic methods. Specific bacterial and parasitic genes were amplified through multiplex PCR and subjected to the NanoCHIP® detection and analysis system. The generated amplicons were electronically addressed to discrete loci on the NanoCHIP® cartridge, pre-activated with specific capture oligonucleotides. Detection was achieved through specific fluorescent reporter oligonucleotides. The molecular method used in the testing laboratory, and STD6 ACE Detection kit of Seegene were used as reference methods.

Results. The NanoCHIP® results were in accordance with the characterizations of the tested samples in terms of clinical sensitivity and specificity. The NanoCHIP® multiplex analysis provided clear results about the identity of the pathogen, either bacterium or parasite, within a working day time frame.

Conclusions. The NanoCHIP® based test has proven to be a useful system for medium-high throughput screening of urine and swab samples for reliably detection of STIs. The NanoCHIP® technology presents significant advantages, mainly in terms of minimal hands-on time, improved laboratory workflow and turnaround time, enabling flexibility and saving costs.

Background

Sexually transmitted infections (STI) are widely abundant all over the world, presenting constantly a yearly growing trend and consequently increased annual health costs. STIs may cause acute symptoms, chronic infections, infertility, ectopic pregnancy, blindness, cervical cancer, susceptibility to HIV infection and death. Most cases are treatable, however no or late treatment may lead to severe health complications. This is further emphasized in view of the common carriage of STIs in asymptomatic state. Therefore, there is a need for detection of early stages of STIs in order to control and prevent these diseases. Savyon Diagnostics has recently developed a novel molecular-based screening test for simultaneous detection of the most abundant STIs utilizing its proprietary NanoCHIP®XL molecular electronic microarray system. The panel of pathogens is composed of *Chlamydia trachomatis*, *Neisseria gonorrhoeae*, *Trichomonas vaginalis*, *Ureaplasma urealyticum/parvum*, *Mycoplasma genitalium* and *Mycoplasma hominis*. The system is able to analyze various specimen types, such as urine, vaginal and urethral swabs, as well as culture isolates. The DNA from these specimens may be extracted by various available automated nucleic acids extraction instruments, but an extraction protocol has been developed in which no instrumentation is needed at the same extraction efficiency. The usage of the NanoCHIP® platform enables high multiplexing capabilities together with testing multiple patient samples in the same run, and thus offers a powerful and unique medium-high throughput screening tool, demonstrating a highly performing, time-saving and cost-effective test.

The test provides:

- Automated testing process from sample to result ; No need for culture enrichment
- Compatibility with the most abundant automatic nucleic acids extraction instruments, as well as with instrument-free extractions
- Versatility: the user may choose to test the complete panel, or any combination of organisms composing the panel
- Detection of bacteria and parasites in the same run
- Minimal hands-on time, vast reduction in the workload of the lab
- Same day results
- Up to 96 samples per run
- Cost-effectiveness

Results

Table 1. STI3 and STI6 Panels Composition

Panel/Pathogen	<i>Chlamydia Trachomatis</i>	<i>Trichomonas Vaginalis</i>	<i>Neisseria Gonorrhoeae</i>	<i>Ureaplasma urealyticum/parvum</i>	<i>Mycoplasma hominis</i>	<i>Mycoplasma genitalium</i>
STI3	+	+	+	-	-	-
STI6	+	+	+	+	+	+

Table 2. Typical STI6 Panel Raw Results and Interpretation

Reference Characterization	Trichomonas vaginalis Reporter	Chlamydia trachomatis Reporter 1	Chlamydia trachomatis Reporter 2	Ureaplasma urealyticum/parvum Reporter	Mycoplasma hominis Reporter	Mycoplasma genitalium Reporter	Neisseria gonorrhoeae Reporter 1	Neisseria gonorrhoeae Reporter 2	Internal Control	NanoCHIP® Interpretation
TV,Uu/Up	45890	622	141	48282	628	823	1640	388	13289	Trichomonas vaginalis, Ureaplasma urealyticum/parvum
TV	8895	589	88	3075	657	602	776	408	5616	Trichomonas vaginalis
TV,Uu/Up	92665	6930	692	37885	237	938	2965	751	13066	Trichomonas vaginalis, Ureaplasma urealyticum/parvum
TV,Uu/Up	68253	926	81	24871	397	315	1938	411	23655	Trichomonas vaginalis, Ureaplasma urealyticum/parvum
Cl	866	74891	21266	1058	601	826	1497	170	14843	Chlamydia trachomatis
Cl	216	64788	14218	2389	427	626	1361	518	14824	Chlamydia trachomatis
Cl,TV,Uu/Up	12065	74199	38699	40220	2084	5763	1523	302	15990	Chlamydia trachomatis, Trichomonas vaginalis, Ureaplasma urealyticum/parvum
Cl,Uu/Up	563	77662	27992	42128	829	837	2242	499	29864	Chlamydia trachomatis, Ureaplasma urealyticum/parvum
Cl,Uu/Up	180	48072	14405	53136	1068	649	2311	443	19995	Chlamydia trachomatis, Ureaplasma urealyticum/parvum
Cl	418	22461	21228	7461	1874	757	1976	375	20377	Chlamydia trachomatis
Uu/Cl	283	13549	141	35886	264	260	1574	1196	15141	Ureaplasma urealyticum/parvum, Chlamydia trachomatis
Uu/Cl	2633	19338	6723	40233	455	842	2352	600	20345	Ureaplasma urealyticum/parvum, Chlamydia trachomatis
Uu/Up	2092	1301	201	17234	2028	584	5803	1657	18954	Ureaplasma urealyticum/parvum
Uu/Up	742	963	470	50357	149	680	1618	287	17495	Ureaplasma urealyticum/parvum
Mh,Cl,Mg	7898	27201	6026	3314	43994	13148	1806	604	9018	Mycoplasma hominis, Chlamydia trachomatis, Mycoplasma genitalium
Mh,TV,Cl	10795	20990	1822	2118	12603	935	1477	439	13258	Mycoplasma hominis, Trichomonas vaginalis, Chlamydia trachomatis
Mh,Uu/Up,Cl	140	20278	9269	18319	51093	885	1633	485	9450	Mycoplasma hominis, Ureaplasma urealyticum/parvum, Chlamydia trachomatis
Mh	2318	673	2492	1364	63883	330	2256	636	13478	Mycoplasma hominis
Mg,Uu/Up	1087	974	197	49308	436	9902	1778	154	19697	Mycoplasma genitalium, Ureaplasma urealyticum/parvum
Mg,Uu/Up	252	4031	109	33055	640	707	26414	8229	13750	Neisseria gonorrhoeae, Ureaplasma urealyticum/parvum
Mg	148	736	213	135	382	1029	30132	301	4656	Neisseria gonorrhoeae
Mg	143	674	182	809	675	1414	6563	16378	19432	Neisseria gonorrhoeae
Mg,Uu/Up	4573	717	109	11032	154	769	53818	17924	18917	Neisseria gonorrhoeae, Ureaplasma urealyticum/parvum
Cutoff	8,000	9,000	9,000	9,000	9,000	7,500	9,000	7,000	5,000	

Samples were characterized according to various RT-PCR based tests that are used in the laboratories that provided the samples
NanoCHIP® results show complete accordance with the reference results

Table 4. STI3 Panel Performance

Pathogen	No. of detected samples / No. of tested samples	% of Detection
<i>Chlamydia trachomatis</i>	50/50	100%
<i>Trichomonas vaginalis</i>	40/40	100%
<i>Neisseria gonorrhoeae</i>	21/21	100%

All the samples were detected in the STI3 system according to their characterization

Table 5. STI6 Panel Performance

Pathogen	No. of detected samples / No. of tested samples	% of Detection
<i>Chlamydia trachomatis</i>	18/18	100%
<i>Trichomonas vaginalis</i>	14/14	100%
<i>Neisseria gonorrhoeae</i>	13/13	100%
<i>Ureaplasma urealyticum/parvum</i>	16/16	100%
<i>Mycoplasma hominis</i>	14/14	100%
<i>Mycoplasma genitalium</i>	12/12	100%

All the samples were detected in the STI6 system according to their characterization

Results (cont.)

Table 3. Typical STI3 Panel Raw Results and Interpretation

Sample #	Reference Characterization	Chlamydia trachomatis Reporter 1	Chlamydia trachomatis Reporter 2	Trichomonas vaginalis Reporter	Neisseria gonorrhoeae Reporter 1	Neisseria gonorrhoeae Reporter 2	Internal control	NanoCHIP Interpretation
1	N	665	237	311	518	986	10,458	Negative
2	N	498	182	263	218	854	12,831	Negative
3	CT	51,644	16,372	226	141	510	10,078	Chlamydia trachomatis
4	CT	67,789	10,967	279	381	563	43,140	Chlamydia trachomatis
5	CT	68,170	29,069	162	475	826	93,389	Chlamydia trachomatis
6	CT	75,431	18,024	205	4,655	646	21,643	Chlamydia trachomatis
7	CT	53,687	41,031	157	390	677	12,517	Chlamydia trachomatis
8	TV	1,464	484	37,384	323	512	23,036	Trichomonas vaginalis
9	TV	2,491	3,059	20,165	359	565	65,746	Trichomonas vaginalis
10	TV	489	2,504	35,627	2,280	504	8,549	Trichomonas vaginalis
11	TV	2,362	346	29,054	885	564	12,293	Trichomonas vaginalis
12	CT	56,079	53,511	284	461	705	7,849	Chlamydia trachomatis
13	CT	65,952	43,989	2,007	244	584	41,340	Chlamydia trachomatis
14	CT	44,917	2,662	1,596	257	571	60,874	Chlamydia trachomatis
15	CT	68,130	8,241	346	1,287	738	43,490	Chlamydia trachomatis
16	TV	1,362	402	49,485	645	762	41,199	Trichomonas vaginalis
17	TV	613	289	54,277	407	897	52,161	Trichomonas vaginalis
18	TV/CT	31,436	5,519	54,264	299	625	65,639	Trichomonas vaginalis, Chlamydia trachomatis
19	TV	822	239	56,515	477	701	76,778	Trichomonas vaginalis
20	TV	541	200	49,258	422	456	56,794	Trichomonas vaginalis
21	TV	1,202	682	54,339	476	589	65,161	Trichomonas vaginalis
22	TV	2,106	476	58,905	499	631	77,459	Trichomonas vaginalis
23	TV	1,068	462	55,322	486	691	76,458	Trichomonas vaginalis
24	TV	760	217	51,152	376	505	63,644	Trichomonas vaginalis
25	TV	943	2,542	55,722	220	580	87,841	Trichomonas vaginalis
26	TV	1,538	314	50,005	874	795	45,914	Trichomonas vaginalis
27	NG	1,765	560	570	60,269	59,972	49,922	Neisseria gonorrhoeae
28	NG	2,236	575	474	59,154	59,768	58,734	Neisseria gonorrhoeae
29	NG	1,486	422	423	634	62,501	18,027	Neisseria gonorrhoeae
30	NG	1,692	383	595	62,853	60,186	19,656	Neisseria gonorrhoeae
31	NG/CT	47,858	12,791	290	20,122	30,975	74,833	Neisseria gonorrhoeae, Chlamydia trachomatis
32	NG	2,018	1,239	242	59,061	56,014	73,292	Neisseria gonorrhoeae
33	NG	504	48	148	57,568	47,948	56,144	Neisseria gonorrhoeae
34	NG	971	190	406	72,805	63,858	43,857	Neisseria gonorrhoeae
35	NG	1,381	609	241	489	59,041	54,490	Neisseria gonorrhoeae
36	NG	1,774	528	246	65,476	70,824	22,166	Neisseria gonorrhoeae
37	NG	606	89	123	50,554	46,786	51,402	Neisseria gonorrhoeae
38	NG	1,026	431	1,054	67,049	71,202	46,808	Neisseria gonorrhoeae
39	NO DNA	1,037	223	375	678	893	818	NO DNA
40	NO DNA	787	191	420	236	563	591	NO DNA
41	NO DNA	882	216	1,744	362	618	1,758	NO DNA
	cutoff	5,000	5,000	10,000	7,000	7,000	3,000	

Samples were characterized according to various RT-PCR based tests that are used in the laboratories that provided the samples

NanoCHIP® results show complete accordance with the reference results

Discussion & Summary

- The current NanoCHIP® system provides two versions of diagnostic panels enabling detection of the most abundant STI pathogens
- The NanoCHIP® results show excellent potency for detection of the target pathogens in both, STI3 as well as STI6 systems, in clinical urethral swabs and in urine specimens, as demonstrated versus various RT-PCR currently used tests
- Both tests are specific to the detected pathogens and there are no cross reactions between the pathogens in the multiplexed assay
- The results are clear and enable definite interpretations to be presented conveniently to the end user
- All in all The NanoCHIP® has proven to be a useful platform for medium-high throughput screening of clinical samples for reliably detection of Sexual Transmitted Infections, and the current tests present significant advantages mainly in terms of minimal hands-on time, improved laboratory workflow and turn around time, enabling flexibility and saving costs

Materials & Methods

DNA was extracted from characterized urethral samples or urine swab specimens using a variety of validated available manual and automatic methods. Specific bacterial and parasitic genes were amplified through multiplex PCR and subjected to the NanoCHIP® detection and analysis system. The generated amplicons were electronically addressed to discrete loci on the NanoCHIP® cartridge, pre-activated with specific capture oligonucleotides. Detection was achieved through specific fluorescent reporter oligonucleotides. The efficiency of the purification was monitored by adding human epithelial cells into each tested sample. Then human β -globin gene was amplified in the same multiplex PCR and analyzed simultaneously in the NanoCHIP® system. The molecular method used in the sample provider laboratory, and STD6 ACE Detection kit of Seegene were used as reference methods

The NanoCHIP® System is an automated multiplex platform capable of detecting and analyzing multiple targets together with multiple samples on the same run utilizing the electronic micro-array technology
More information on the system and its function can be found at Booth #81 and at:
www.nanochipxl.com and www.savyondx.com



Objective

The aim of this work is to present the newly developed tests and their utility for screening purposes of Sexual Transmitted Infections especially in high-risk populations

The NanoCHIP® Workflow

