

eP321

ePoster Viewing

Diagnostics: detection of ESBLs and carbapenemases

PILOT EVALUATION OF THE BLACTA TEST FOR PREDICTING 3RD GENERATION CEPHALOSPORIN RESISTANCE (3GCR) IN ESCHERICHIA COLI, KLEBSIELLA PNEUMONIAE, KLEBSIELLA OXYTOCA AND PROTEUS MIRABILIS FROM BLOOD

B.M. Willey¹, X. Trimi¹, P. Lo¹, S.M. Poutanen²

¹Microbiology, Mount Sinai Hospital/University Health Network, Toronto Ontario, Canada ; ²Microbiology, Mount Sinai Hospital/University Health Network/University of Toronto, Toronto Ontario, Canada

Objectives: VITEK MS (VMS; bioMérieux) testing of blood culture isolates subcultured and incubated for ~3-6h on agar (short-incubation blood culture isolates) (SIBCI) significantly improved time-to-reporting of organism identifications (ID). No similar improvement in antimicrobial susceptibility test (AST) reporting has been achieved. We conducted a pilot evaluation of the BLACTA test (BLT; Bio-Rad) for predicting 3GCR in SIBCI *E. coli* (EC), *K. pneumoniae* (KP), *K. oxytoca* (KO), and *P. mirabilis* (PM). The BLT is a low complexity method that detects 3GCR from isolates within 15-minutes by enzymatic hydrolysis of a chromogenic cephalosporin resulting in a yellow to red colour-change.

Methods: 153 prospective semi-consecutive SIBCI were tested by BLT and interpreted as per kit instructions immediately following VMS ID as EC, KP, KO or PM. BLT results were compared 24-48h later to VITEK 2 N213 AST, Double Disk Test results, and molecular detection of b-lactamases if available. BLT-sensitivities and specificities for detecting class A extended spectrum b-lactamases (A-ESBL) and for 3GCR overall were calculated and the impact on 3GCR reporting-time was compared to conventional AST.

Results: Of 153 BLT completed, 103 (66.9%) were completed on EC, 37 (24%) on KP, 5 (3.3%) on KO, and 8 (5.2%) on PR. Of these, 22 (14.3%) were BLT-positive (red), 1 (0.7%) was indeterminate (orange), and 131 (85.1%) were negative (yellow). Of the 22 BLT-positives, 100% (95%CI: 84.2-100) were 3GCR; 11 were A-ESBL EC (7 CTX-M, 4 TEM/SHV), 7 were A-ESBL KP (4 CTX-M, 3 TEM/SHV), 2 were A-ESBL KO (2 CTX-M), while 2 KO were hyper-producers of chromosomal class A K1 (OXY) b-lactamase. The single BLT-indeterminate and 4 of the 131 (3.1%) BLT-negative isolates were 3GCR due to Class C b-lactamases; all 5 were multidrug-R EC. Sensitivities/specificities for class A-ESBL and for 3GCR overall were 100% (95%CI: 82.5-100)/98.5% (94.4-99.9) and 81.5% (95%CI: 62.8-92.3)/100% (96.4-100), respectively. Compared to the conventional AST, 3GCR was reported based on BLT-positive results 24h earlier in 83.3% of cases (P<0.0001). In addition, BLT identified an A-ESBL (TEM/SHV) KP that would have been missed on conventional AST due to a mixed culture with overgrowth of a 3CG susceptible KP.

Conclusions: In this pilot study, the BLT was found to be rapid and low-complexity method that proved to be highly accurate in its prediction of A-ESBL and 3GCR from SIBCI. Through interim reporting of BLT-positive results, patients would have had the benefit of appropriate treatment for 3GCR organisms 24h earlier than in the current conventional AST algorithm, an improvement considered extremely significant.