

**P0950**

**Paper Poster Session**

**Recent studies of non-culture techniques for detection of resistance**

**Rapid carbapenemase- and ESBL-producing strains detection from positive blood culture by using CarbaNP and ESBL-NDP tests**

Elena Pilotto<sup>1</sup>, Laura Maccacaro<sup>2</sup>, Manuela Pegoraro<sup>2</sup>, Giuseppe Cornaglia<sup>3</sup>, Annarita Mazzariol<sup>\*3</sup>

<sup>1</sup>*Microbiology, Diagnostics and Public Health, Verona, Italy*

<sup>2</sup>*Microbiology and Virology, Azienda Ospedaliera Universitaria Integrata Verona, Pathology and Diagnostics, Verona, Italy*

<sup>3</sup>*University of Verona, Diagnostics and Public Health, Verona, Italy*

**Background:** This study aims at optimizing the application of CarbaNP test and ESBL NDP test directly from positive blood culture samples, with a new flowchart that reduces to at few hours the turnaround time for identification of carbapenemase- and ESBL-producing strains.

**Material/methods:** A total of 35 positive blood cultures were included in the study, namely 22 with ESBL-producing strains, and 13 with carbapenemase-producing strains. The study included *Klebsiella pneumoniae*, *Escherichia coli*, and *Pseudomonas aeruginosa*.

All isolates underwent identification by Maldi-TOF, antimicrobial susceptibility test by microdilution, PCR for detection of ESBL and carbapenemases genes, carbaNP and ESBL-NDP tests.

CarbaNP and ESBL-NDP were applied to positive blood cultures under different sub-cultures conditions, namely: amount of blood to be subcultured (3 drops vs 10 drops); type of culture media (Muller-Hinton agar vs blood agar); incubation conditions (air vs 5% of CO<sub>2</sub>). All subsets were evaluated for the CarbaNP and ESBL-NDP positivity at 2, 3 and 4 hours of incubation.

**Results:** 10 blood drops of positive blood cultures harbouring carbapenemases producing strains streaked on MH plates and incubated at 37°C on air yielded a 61% positivity of carbaNP test after 2 h, as compared to 23% with 3 drops; in both cases positivity raised up to 92% when the test was performed after 4 hours.

3 blood drops incubated at 37°C on air yielded after 2 hours a 61% and 23% positive CarbaNP test for blood agar and MH, respectively. At 4 hours positivity was 100% for blood agar and 92% for MH.

The evaluation of air incubation vs air with 5% of CO<sub>2</sub>, performed on blood agar with 10 drops (i.e., the best conditions stemming from the previous results) showed a 100% of positive carbaNP test at 3 hours for air plus 5% of CO<sub>2</sub>, compared to 92% only on air, and a 100% positivity at 4 hours for both

Results for all beta-lactamase producing strains (tested with the ESBL NDP test) paralleled those described above for the strains resistant to carbapenems.

**Conclusions:** The Carba NP and ESBL NDP test were optimized for direct performance from positive blood cultures.

The best results were obtained with 10 drops of positive blood culture streaked on blood agar plates and incubated at 37°C for up to 3 hours in air added with 5% of CO<sub>2</sub>. The carbaNP and ESBL-NDP test were applied directly to positive blood culture in a lab routine setting, and reduced the reporting time to a 5 hours.

Noteworthy, the same subculture can be used for Maldi-TOF identification and antimicrobial susceptibility test, too, thus allowing one to obtain a complete report in 24 h.